

CAB21164

Product Information

Product SKU:	CAB21164	Gene ID:	161176	Size:	20uL, 100uL
Clone No:	ARC3013	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	112kDa	Conjugate:	Unconjugated
Calculated MW:	112kDa	Isotype:	IgG

Immunogen Information

Background: Enables actin filament binding activity and cytoskeleton-nuclear membrane anchor activity. Involved in cytoskeleton organization; establishment of protein localization to membrane; and regulation of cell shape. Located in nuclear membrane. Part of meiotic nuclear membrane microtubule tethering complex. Biomarker of Huntington's disease.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: KASH3; NET53; Nesp3; C14orf49; C14orf139; LINC00341; NCRNA00341; Nesprin3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 750-935 of human Nesprin3 (NP_689805.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.