

CAB21189

Product Information

Product SKU:	CAB21189	Gene ID:	6856	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	29kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in chemical synaptic transmission. Located in extracellular exosome.

Recommended Dilution: IF/ICC, 1:50 - 1:200

Synonyms: SYPL; H-SP1; SYPL1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 50-120 of human SYPL1 (NP_006745.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.