

CAB21210

Product Information

Product SKU:	CAB21210	Gene ID:	9868	Size:	20uL, 100uL
Clone No:	ARC53258	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	75kDa	Conjugate:	Unconjugated
Calculated MW:	67kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes an import receptor of the outer mitochondrial membrane that is part of the translocase of the outer membrane complex. This protein is involved in the import of mitochondrial precursor proteins. The protein interacts with the SARS-CoV and SARS-Cov-2 ORF9b proteins.

Recommended Dilution: WB,1:1000 - 1:5000

Synonyms: Tom70; TOMM70A; TOM70

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 245-608 of human TOM70 (NP_055635.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.