

## CAB21214

---

### Product Information

<b>Product SKU:</b>	CAB21214	<b>Gene ID:</b>	8290/8350	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC53046	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat,Other (Wide Range Predicted)

---

### Additional Information

<b>Observed MW:</b>	17kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	15kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. This structure consists of approximately 146 bp of DNA wrapped around a nucleosome, an octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H3 family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. This gene is found in the large histone gene cluster on chromosome 6p22-p21.3.

**Recommended Dilution:** WB,1:1000 - 1:5000

**Synonyms:** H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Lactic acid-Histone H3-K18

**Purification Method:** Affinity purification

**Immunogen:** A synthetic lactylated peptide around K18 of human Lactic acid-Histone H3-K18 (NP\_003520.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.