

## CAB21244

---

### Product Information

<b>Product SKU:</b>	CAB21244	<b>Gene ID:</b>	7040	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC5123-07	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	45-60kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	44kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate a latency-associated peptide (LAP) and a mature peptide, and is found in either a latent form composed of a mature peptide homodimer, a LAP homodimer, and a latent TGF-beta binding protein, or in an active form consisting solely of the mature peptide homodimer. The mature peptide may also form heterodimers with other TGFbeta family members. This encoded protein regulates cell proliferation, differentiation and growth, and can modulate expression and activation of other growth factors including interferon gamma and tumor necrosis factor alpha. This gene is frequently upregulated in tumor cells, and mutations in this gene result in Camurati-Engelmann disease.

**Recommended Dilution:** WB,1:2000 - 1:10000

**Synonyms:** CED; LAP; DPD1; TGFB; IBDIMDE; TGFbeta; TGF-beta1; TGF beta 1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of human TGF beta 1.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.