

CAB21255

Product Information

Product SKU:	CAB21255	Gene ID:	9451	Size:	20uL, 100uL
Clone No:	ARC52369	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	140kDa	Conjugate:	Unconjugated
Calculated MW:	125kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene phosphorylates the alpha subunit of eukaryotic translation-initiation factor 2, leading to its inactivation, and thus to a rapid reduction of translational initiation and repression of global protein synthesis. This protein is thought to modulate mitochondrial function. It is a type I membrane protein located in the endoplasmic reticulum (ER), where it is induced by ER stress caused by malformed proteins. Mutations in this gene are associated with Wolcott-Rallison syndrome.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: PEK; WRS; PERK

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 580-720 of human PERK (NP_004827.4).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.