

CAB21273

Product Information

Product SKU:	CAB21273	Gene ID:	838804	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Arabidopsis thaliana

Additional Information

Observed MW:	88kDa	Conjugate:	Unconjugated
Calculated MW:	88kDa	Isotype:	IgG

Immunogen Information

Background: encodes a serine/threonine protein phosphatase with an N-terminal Kelch-repeat domain, which is nuclear localized and expressed preferentially in elongating cells. Genetic evidence suggest that this gene plays a redundant role (along with other members of the same gene family) in modulating growth in response to brassinosteroid.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: BRI1 SUPPRESSOR 1; F21B7.7; BSU1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 694-793 of arabidopsis thaliana BSU1 (NP_171844.6).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.