

CAB21292

Product Information

Product SKU:	CAB21292	Gene ID:	54918	Size:	20uL, 100uL
Clone No:	ARC53656	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	20-22kDa	Conjugate:	Unconjugated
Calculated MW:	20kDa	Isotype:	IgG

Immunogen Information

Background: This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and transmembrane 4 superfamilies. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 3. This gene is widely expressed in many tissues, but the exact function of the encoded protein is unknown.

Recommended Dilution: WB,1:2000 - 1:20000

Synonyms: CKLFSF6; PRO2219; CMTM6

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 84-183 of human CMTM6 (NP_060271.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.