

CAB21294

Product Information

Product SKU:	CAB21294	Gene ID:	9276	Size:	20uL, 100uL
Clone No:	ARC52626	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	102kDa	Conjugate:	Unconjugated
Calculated MW:	102kDa	Isotype:	IgG

Immunogen Information

Background: The Golgi coatomer complex (see MIM 601924) constitutes the coat of nonclathrin-coated vesicles and is essential for Golgi budding and vesicular trafficking. It consists of 7 protein subunits, including COPB2.

Recommended Dilution: WB,1:3000 - 1:20000

Synonyms: OPDD; MCPH19; beta'-COP; COPB2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 657-906 of human COPB2 (NP_004757.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.