

## CAB21297

---

### Product Information

<b>Product SKU:</b> CAB21297	<b>Gene ID:</b> 5452	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b> 60-75kDa	<b>Conjugate:</b> -
<b>Calculated MW:</b> 51kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** The protein encoded by this gene is a homeobox-containing transcription factor of the POU domain family. The encoded protein binds the octamer sequence 5'-ATTGTCAT-3', a common transcription factor binding site in immunoglobulin gene promoters. Several transcript variants encoding different isoforms have been found for this gene.

**Recommended Dilution:** WB,1:100 - 1:500

**Synonyms:** OCT2; OTF2; Oct-2; POU2F2

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence within amino acids 400-479 of human POU2F2 (NP\_001193954.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.