

CAB21303

Product Information

Product SKU:	CAB21303	Gene ID:	79037	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	30-40kd	Conjugate:	-
Calculated MW:	34kDa	Isotype:	IgG

Immunogen Information

Background: Enables phosphatase binding activity and signaling receptor activity. Involved in negative regulation of T cell receptor signaling pathway. Located in plasma membrane.

Recommended Dilution: WB, 1:100 - 1:500

Synonyms: CD112R; C7orf15; PVRIG

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 193-326 of human PVRIG (NP_076975.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.