

## CAB21350

---

### Product Information

<b>Product SKU:</b>	CAB21350	<b>Gene ID:</b>	11168	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	Refer to figures	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	60kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Enables DNA-binding transcription factor binding activity; chromatin binding activity; and transcription coactivator activity. Involved in mRNA 5'-splice site recognition and positive regulation of transcription by RNA polymerase II. Located in heterochromatin; nuclear periphery; and nucleoplasm.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** p52; p75; PAIP; DFS70; LEDGF; PSIP2; PSIP1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 301-530 of human PSIP1 (NP\_150091.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.