

## CAB21351

---

### Product Information

<b>Product SKU:</b> CAB21351	<b>Gene ID:</b> 11171	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b> 39kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 38kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** Enables RNA binding activity. Involved in maintenance of gastrointestinal epithelium; negative regulation of transforming growth factor beta receptor signaling pathway; and spliceosomal snRNP assembly. Located in cytosol. Part of SMN complex. Implicated in adenocarcinoma; colorectal carcinoma; large cell carcinoma; lung carcinoma; and squamous cell neoplasm.

**Recommended Dilution:** WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200 IP,0.5µg-4µg antibody for 200µg-400µg extracts of whole cells

**Synonyms:** MAWD; PT-WD; UNRIP; AP

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human STRAP (NP\_009109.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.