

CAB21352

Product Information

Product SKU:	CAB21352	Gene ID:	11252	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	56kDa	Isotype:	IgG

Immunogen Information

Background: This gene is a member of the protein kinase C and casein kinase substrate in neurons family. The encoded protein is involved in linking the actin cytoskeleton with vesicle formation by regulating tubulin polymerization. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: SDPII; PACSIN2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 240-410 of human PACSIN2 (NP_009160.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.