

CAB21388

Product Information

Product SKU:	CAB21388	Gene ID:	200081	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	75kDa	Conjugate:	Unconjugated
Calculated MW:	62kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable syntaxin binding activity. Predicted to be involved in exocytosis. Predicted to act upstream of or within B cell activation. Located in cytoplasm.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: IL14; TXLN; TXLNA

Purification Method: Affinity purification

Immunogen: Recombinant Protein corresponding to a sequence within amino acids 1-160 of human TXLNA (NP_787048.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.