

CAB21407

Product Information

Product SKU: CAB21407	Gene ID: 23170	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: Refer to figures	Conjugate: Unconjugated
Calculated MW: 74kDa	Isotype: IgG

Immunogen Information

Background: Enables H4K20me3 modified histone binding activity and tubulin binding activity. Involved in negative regulation of type I interferon-mediated signaling pathway and regulation of mitotic cell cycle. Located in cytoplasm.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: dJ526I14.2; TTLL12

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 485-644 of human TTLL12 (NP_055955.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.