

CAB21408

Product Information

Product SKU:	CAB21408	Gene ID:	2319	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	47kDa	Conjugate:	Unconjugated
Calculated MW:	47kDa	Isotype:	IgG

Immunogen Information

Background: Caveolae are small domains on the inner cell membrane involved in vesicular trafficking and signal transduction. This gene encodes a caveolae-associated, integral membrane protein, which is thought to function in neuronal signaling.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: ESA; ECS1; ESA1; ECS-1; M17S1; FLOT2

Purification Method: Affinity purification

Immunogen: Recombinant Protein corresponding to a sequence within amino acids 239-428 of human FLOT2 (NP_004466.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.