

CAB21432

---

## Product Information

<b>Product SKU:</b>	CAB21432	<b>Gene ID:</b>	27037	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

## Additional Information

<b>Observed MW:</b>	Refer to figures	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	69kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** The protein encoded by this gene is of unknown function. However, it is orthologous to the mouse Trmt2a gene and contains an RNA methyltransferase domain. Expression of this gene varies during the cell cycle, with aberrant expression being a possible biomarker in certain breast cancers. Several transcript variants encoding two different isoforms have been found for this gene.

**Recommended Dilution:** WB, 1:500 - 1:2000 IHC-P, 1:50 - 1:200

**Synonyms:** HTF9C; TRMT2A/HTF9C

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human TRMT2A/HTF9C (NP\_073564.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.