

CAB21435

Product Information

Product SKU:	CAB21435	Gene ID:	27250	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	58kDa	Conjugate:	Unconjugated
Calculated MW:	52kDa	Isotype:	IgG

Immunogen Information

Background: This gene is a tumor suppressor and encodes a protein that binds to the eukaryotic translation initiation factor 4A1 and inhibits its function by preventing RNA binding. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: H731; PDCD4

Purification Method: Affinity purification

Immunogen: Recombinant Protein corresponding to a sequence within amino acids 1-260 of human PDCD4 (NP_055271.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.