

CAB21438

Product Information

Product SKU: CAB21438	Gene ID: 2810	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: 27kDa	Conjugate: Unconjugated
Calculated MW: 28kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a cell cycle checkpoint protein. The encoded protein binds to translation and initiation factors and functions as a regulator of mitotic translation. In response to DNA damage this protein plays a role in preventing DNA errors during mitosis.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: YWHAS; ma

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-248 of human 14-3-3 sigma (NP_006133.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.