

CAB21439

Product Information

Product SKU:	CAB21439	Gene ID:	283373	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	115kDa	Isotype:	IgG

Immunogen Information

Background: Putative regulatory subunit of protein phosphatase 6 (PP6 that may be involved in the recognition of phosphoprotein substrates.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: ANKRD33; ANKRD52

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 450-660 of human ANKRD52 (NP_775866.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.