

## CAB21481

---

### Product Information

<b>Product SKU:</b>	CAB21481	<b>Gene ID:</b>	3916	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	90-120kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	45kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may also play a role in tumor cell metastasis.

**Recommended Dilution:** WB, 1:500 - 1:1000 IF/ICC, 1:50 - 1:200

**Synonyms:** LAMPA; CD107a; LGP120; P1

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 29-382 of human LAMP1 (NP\_005552.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.