

## CAB21483

---

### Product Information

<b>Product SKU:</b>	CAB21483	<b>Gene ID:</b>	3956	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	14kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	15kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. This gene product may act as an autocrine negative growth factor that regulates cell proliferation.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** GBP; GAL1; Galectin 1/LGALS1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant Protein corresponding to a sequence within amino acids 1-135 of human Galectin 1/LGALS1 (NP\_002296.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.