

## CAB21498

---

### Product Information

<b>Product SKU:</b>	CAB21498	<b>Gene ID:</b>	4580	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

---

### Additional Information

<b>Observed MW:</b>	Refer to figures	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	51kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Predicted to be involved in mitochondrion organization. Part of MIB complex and SAM complex.

**Recommended Dilution:** WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

**Synonyms:** MTX; MTXN; MTX1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 186-435 of human MTX1 (NP\_942584.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.