

CAB21523

Product Information

Product SKU: CAB21523	Gene ID: 50674	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Mouse

Additional Information

Observed MW: Refer to figures	Conjugate: Unconjugated
Calculated MW: 23kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene is a basic helix-loop-helix (bHLH) transcription factor involved in neurogenesis. The encoded protein likely acts as a heterodimer with another bHLH protein. Defects in this gene are a cause of congenital malabsorptive diarrhea 4 (DIAR4).

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: ngn3; Atoh5; NGN-3; Math4B; bHLHa7; NEUROG3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 144-214 of human NEUROG3 (NP_066279.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.