

## CAB21530

---

### Product Information

<b>Product SKU:</b>	CAB21530	<b>Gene ID:</b>	51070	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	33kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	33kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** The protein encoded by this gene may modulate the activity and localization of nitric oxide synthase (endothelial and neuronal) and thus nitric oxide production. Alternative splicing results in multiple transcript variants that encode the same protein.

**Recommended Dilution:** WB, 1:100 - 1:500

**Synonyms:** CGI-25; NOSIP

**Purification Method:** Affinity purification

**Immunogen:** Recombinant Protein corresponding to a sequence within amino acids 1-301 of human NOSIP (NP\_057037.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.