

CAB21538

Product Information

Product SKU: CAB21538	Gene ID: 5214	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: Refer to figures	Conjugate: Unconjugated
Calculated MW: 86kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a member of the phosphofructokinase A protein family. The encoded enzyme is the platelet-specific isoform of phosphofructokinase and plays a key role in glycolysis regulation. This gene may play a role in metabolic reprogramming in some cancers, including clear cell renal cell carcinomas, and cancer of the bladder, breast, and lung. Alternative splicing results in multiple transcript variants.

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	PFKE; PFK-C; PFK-P; ATP-PFK; PFKP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 605-784 of human PFKP (NP_002618.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.