

## CAB21539

---

### Product Information

<b>Product SKU:</b> CAB21539	<b>Gene ID:</b> 5217	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b> 15kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 15kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** The protein encoded by this gene is a ubiquitous actin monomer-binding protein belonging to the profilin family. It is thought to regulate actin polymerization in response to extracellular signals. There are two alternatively spliced transcript variants encoding different isoforms described for this gene.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** PFL; D3S1319E; PFN2

**Purification Method:** Affinity purification

**Immunogen:** Recombinant Protein corresponding to a sequence within amino acids 1-140 of human PFN2 (NP\_444252.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.