

CAB21555

Product Information

Product SKU:	CAB21555	Gene ID:	5481	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	39kDa	Conjugate:	Unconjugated
Calculated MW:	41kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a member of the peptidyl-prolyl cis-trans isomerase (PPIase) family. PPIases catalyze the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and accelerate the folding of proteins. This protein has been shown to possess PPIase activity and, similar to other family members, can bind to the immunosuppressant cyclosporin A.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: CYPD; CYP-40; Cyclophilin 40

Purification Method: Affinity purification

Immunogen: Recombinant Protein corresponding to a sequence within amino acids 161-370 of human Cyclophilin 40 (NP_005029.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.