

CAB21572

Product Information

Product SKU: CAB21572	Gene ID: 56339	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 70kDa	Conjugate: Unconjugated
Calculated MW: 64kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes the 70 kDa subunit of MT-A which is part of N6-adenosine-methyltransferase. This enzyme is involved in the posttranscriptional methylation of internal adenosine residues in eukaryotic mRNAs, forming N6-methyladenosine.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:10 - 1:20 IHC-P,1:500 - 1:1000

Synonyms: M6A; IME4; Spo8; MT-A70; hMETTL3; METTL3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of human METTL3 (NP_062826.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.