

CAB21584

Product Information

Product SKU: CAB21584	Gene ID: 5879	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: Refer to figures	Conjugate: Unconjugated
Calculated MW: 21kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene is a GTPase which belongs to the RAS superfamily of small GTP-binding proteins. Members of this superfamily appear to regulate a diverse array of cellular events, including the control of cell growth, cytoskeletal reorganization, and the activation of protein kinases. Two transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: IF/ICC, 1:50 - 1:200

Synonyms: MIG5; MRD48; Rac-1; TC-25; p21-Rac1; Rac1/2/3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-192 of human Rac1/2/3 (NP_008839.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.