

CAB21659

Product Information

Product SKU:	CAB21659	Gene ID:	79705	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	250kDa	Conjugate:	Unconjugated
Calculated MW:	225kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a multi-domain protein that is a leucine-rich repeat kinase and a GDP/GTP binding protein. The encoded protein is thought to play a role in the regulation of bone mass. Mice lacking a similar gene showed severe osteopetrosis, increased bone mineralization and decreased bone resorption.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: OSMD; RIPK6; Roco1; LRRK1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 10-90 of human LRRK1 (NP_078928.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.