

CAB21679

Product Information

Product SKU:	CAB21679	Gene ID:	8408	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	140kDa	Conjugate:	Unconjugated
Calculated MW:	113kDa	Isotype:	IgG

Immunogen Information

Background: Enables identical protein binding activity; protein serine/threonine kinase activity; and small GTPase binding activity. Involved in several processes, including autophagosome assembly; positive regulation by symbiont of host autophagy; and protein phosphorylation. Located in autophagosome; cytosol; and phagophore assembly site membrane. Is extrinsic component of autophagosome membrane; extrinsic component of omegasome membrane; and extrinsic component of phagophore assembly site membrane. Part of Atg1/ULK1 kinase complex.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: ATG1; ATG1A; UNC51; hATG1; Unc51.1; ULK1

Purification Method: Affinity purification

Immunogen: Recombinant protein of human ULK1.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.