

CAB21687

Product Information

Product SKU: CAB21687	Gene ID: 8454	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 90kDa	Conjugate: Unconjugated
Calculated MW: 90kDa	Isotype: IgG

Immunogen Information

Background: Predicted to enable ubiquitin protein ligase binding activity and ubiquitin-protein transferase activity. Involved in SCF-dependent proteasomal ubiquitin-dependent protein catabolic process and protein ubiquitination. Located in plasma membrane. Part of Parkin-FBXW7-Cul1 ubiquitin ligase complex and SCF ubiquitin ligase complex.

Recommended Dilution: WB,1:100 - 1:500 IHC-P,1:50 - 1:200

Synonyms: CUL1; MGC149834; MGC149835; cullin-1; Cullin 1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 435-593 of human Cullin 1 (NP_003583.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.