

CAB21711

Product Information

Product SKU: CAB21711	Gene ID: 9352	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: Refer to figures	Conjugate: Unconjugated
Calculated MW: 32kDa	Isotype: IgG

Immunogen Information

Background: Enables disulfide oxidoreductase activity. Located in cytosol.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: TxI; TXNL; TRP32; TXL-1; HEL-S-114; TXNL1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-289 of human TXNL1 (NP_004777.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.