

## CAB21745

---

### Product Information

<b>Product SKU:</b>	CAB21745	<b>Gene ID:</b>	3488	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	31kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	31kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Enables insulin-like growth factor I binding activity. Involved in several processes, including cellular response to cAMP; regulation of smooth muscle cell migration; and regulation of smooth muscle cell proliferation. Part of insulin-like growth factor ternary complex. Biomarker of pulmonary fibrosis.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** IBP5; P5

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 21-272 of human IGFBP5 (NP\_000590.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.