

CAB21751

Product Information

Product SKU: CAB21751	Gene ID: 4846	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: Refer to figures	Conjugate: Unconjugated
Calculated MW: 133kDa	Isotype: IgG

Immunogen Information

Background: Nitric oxide is a reactive free radical which acts as a biologic mediator in several processes, including neurotransmission and antimicrobial and antitumoral activities. Nitric oxide is synthesized from L-arginine by nitric oxide synthases. Variations in this gene are associated with susceptibility to coronary spasm. Alternative splicing and the use of alternative promoters results in multiple transcript variants.

Recommended Dilution: IF/ICC, 1:50 - 1:200

Synonyms: eNOS; ECNOS

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1035-1110 of human eNOS (NP_000594.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.