

CAB21760

Product Information

Product SKU:	CAB21760	Gene ID:	84962	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	55kDa	Conjugate:	Unconjugated
Calculated MW:	57kDa	Isotype:	IgG

Immunogen Information

Background: Enables alpha-catenin binding activity and transcription corepressor activity. Involved in several processes, including negative regulation of hippo signaling; positive regulation of gene silencing by miRNA; and regulation of cellular response to hypoxia. Acts upstream of or within gene silencing by miRNA and positive regulation of protein-containing complex assembly. Located in several cellular components, including Golgi apparatus; P-body; and nucleoplasm.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: JUB; ba

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 201-290 of human Ajuba (NP_116265.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.