

CAB21794

Product Information

Product SKU:	CAB21794	Gene ID:	6935	Size:	20uL, 100uL
Clone No:	ARC53599	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	200kDa	Conjugate:	Unconjugated
Calculated MW:	124kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a zinc finger transcription factor. The encoded protein likely plays a role in transcriptional repression of interleukin 2. Mutations in this gene have been associated with posterior polymorphous corneal dystrophy-3 and late-onset Fuchs endothelial corneal dystrophy. Alternatively spliced transcript variants encoding different isoforms have been described.

Recommended Dilution: WB,1:1000 - 1:5000

Synonyms: BZP; TCF8; AREB6; FECD6; NIL2A; PPCD3; ZFHEP; ZFH1A; DELTAEF1; ZEB1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 800-900 of human ZEB1 (NP_110378.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.