

CAB21799

Product Information

Product SKU:	CAB21799	Gene ID:	3326	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	83kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the heat shock protein 90 family; these proteins are involved in signal transduction, protein folding and degradation and morphological evolution. This gene encodes the constitutive form of the cytosolic 90 kDa heat-shock protein and is thought to play a role in gastric apoptosis and inflammation. Alternative splicing results in multiple transcript variants. Pseudogenes have been identified on multiple chromosomes.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: HSP84; HSPC2; HSPCB; D6S182; HSP90B; Hsp90 β

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human Hsp90 β (NP_031381.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.09% Sodium azide,50% glycerol,pH7.3.