

## CAB21809

---

### Product Information

<b>Product SKU:</b>	CAB21809	<b>Gene ID:</b>	2625	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse

---

### Additional Information

<b>Observed MW:</b>	48kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	48kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a protein which belongs to the GATA family of transcription factors. The protein contains two GATA-type zinc fingers and is an important regulator of T-cell development and plays an important role in endothelial cell biology. Defects in this gene are the cause of hypoparathyroidism with sensorineural deafness and renal dysplasia.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** HDR; HDRS; GATA3

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human GATA3 (NP\_001002295.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.