

CAB21821

Product Information

Product SKU: CAB21821	Gene ID: 1069	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 22kDa	Conjugate: Unconjugated
Calculated MW: 20kDa	Isotype: IgG

Immunogen Information

Background: Caltractin belongs to a family of calcium-binding proteins and is a structural component of the centrosome. The high level of conservation from algae to humans and its association with the centrosome suggested that caltractin plays a fundamental role in the structure and function of the microtubule-organizing center, possibly required for the proper duplication and segregation of the centrosome.

Recommended Dilution: WB,1:1000 - 1:5000 IF/ICC,1:50 - 1:200

Synonyms: CALT; CEN2; CETN2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human CETN2 (NP_004335.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.