

CAB21824

---

## Product Information

<b>Product SKU:</b>	CAB21824	<b>Gene ID:</b>	1889	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	125-135kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	87kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** The protein encoded by this gene is involved in proteolytic processing of endothelin precursors to biologically active peptides. Mutations in this gene are associated with Hirschsprung disease, cardiac defects and autonomic dysfunction. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

**Recommended Dilution:** WB,1:500 - 1:2000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

**Synonyms:** ECE; ECE1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 670-770 of human ECE1 (NP\_001388.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.