

CAB21828

Product Information

Product SKU: CAB21828	Gene ID: 3275	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: 49kDa	Conjugate: Unconjugated
Calculated MW: 49kDa	Isotype: IgG

Immunogen Information

Background: Enables several functions, including nuclear receptor binding activity; peroxisome proliferator activated receptor binding activity; and protein homodimerization activity. Involved in several processes, including histone methylation; regulation of androgen receptor signaling pathway; and regulation of transcription, DNA-templated. Located in cytosol and nucleoplasm.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: HRMT1L1; PRMT2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 290-433 of human PRMT2 (NP_001526.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.