

CAB21833

---

## Product Information

<b>Product SKU:</b>	CAB21833	<b>Gene ID:</b>	8407	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

---

## Additional Information

<b>Observed MW:</b>	Refer to figures	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	22kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** The protein encoded by this gene is similar to the protein transgelin, which is one of the earliest markers of differentiated smooth muscle. The specific function of this protein has not yet been determined, although it is thought to be a tumor suppressor. Multiple transcript variants encoding different isoforms have been found for this gene.

**Recommended Dilution:** WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200

**Synonyms:** HA1756; TAGLN2

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-158 of human TAGLN2 (NP\_003555.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.