

CAB21927

Product Information

Product SKU:	CAB21927	Gene ID:	7942	Size:	20uL, 100uL
Clone No:	ARC53778	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	53kDa	Isotype:	IgG

Immunogen Information

Background: Enables DNA-binding transcription factor activity; enzyme binding activity; and transcription cis-regulatory region binding activity. Involved in several processes, including cellular response to amino acid starvation; lysosome localization; and positive regulation of autophagy. Located in cytosol; lysosomal membrane; and nucleoplasm.

Recommended Dilution: IHC-P,1:50 - 1:200

Synonyms: TCFEB; BHLHE35; ALPHATFEB; TFEB

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 50-173 of human TFEB (NP_009093.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.