

CAB21949

Product Information

Product SKU: CAB21949	Gene ID: 9076	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Mouse,Rat

Additional Information

Observed MW: 23kDa	Conjugate: Unconjugated
Calculated MW: 23kDa	Isotype: IgG

Immunogen Information

Background: Tight junctions represent one mode of cell-to-cell adhesion in epithelial or endothelial cell sheets, forming continuous seals around cells and serving as a physical barrier to prevent solutes and water from passing freely through the paracellular space. These junctions are comprised of sets of continuous networking strands in the outwardly facing cytoplasmic leaflet, with complementary grooves in the inwardly facing extracytoplasmic leaflet. The protein encoded by this gene, a member of the claudin family, is an integral membrane protein and a component of tight junction strands. Loss of function mutations result in neonatal ichthyosis-sclerosing cholangitis syndrome.

Recommended Dilution: WB,1:100 - 1:500

Synonyms: CLD1; SEMP1; ILVASC; Claudin 1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 137-211 of human Claudin 1 (NP_066924.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.