

## CAB21967

---

### Product Information

<b>Product SKU:</b>	CAB21967	<b>Gene ID:</b>	10495	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC53634	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Rat

---

### Additional Information

<b>Observed MW:</b>	65-70kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	70kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene is a tumor-specific member of the ECTO-NOX family of genes that encode cell surface NADH oxidases. The encoded protein has two enzymatic activities: catalysis of hydroquinone or NADH oxidation, and protein disulfide interchange. The protein also displays prion-like properties. Alternative splicing results in multiple transcript variants.

**Recommended Dilution:** WB,1:1000 - 1:5000

**Synonyms:** APK1; tNOX; COVA1; ENOX2

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 411-610 of human ENOX2 (NP\_872114.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.