

CAB21994

Product Information

Product SKU: CAB21994	Gene ID: 3297	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse

Additional Information

Observed MW: 73kDa	Conjugate: Unconjugated
Calculated MW: 52kDa/57kDa	Isotype: IgG

Immunogen Information

Background: The product of this gene is a heat-shock transcription factor. Transcription of heat-shock genes is rapidly induced after temperature stress. Hsp90, by itself and/or associated with multichaperone complexes, is a major repressor of this gene.

Recommended Dilution: WB,1:50 - 1:200 IHC-P,1:50 - 1:200

Synonyms: HSF1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 367-529 of human HSF1 (NP_005517.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.