

CAB22001

Product Information

Product SKU:	CAB22001	Gene ID:	139341	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	17kDa	Conjugate:	Unconjugated
Calculated MW:	17kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a protein with a FUN14 superfamily domain. The function of the encoded protein is not known.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: FUNDC1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 50-155 of human FUNDC1 (NP_776155.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.